## **TANH Equation**

The TANH equation describes the relationship between the partial pressure of oxygen (p) and saturation of haemoglobin (s) in human blood.<sup>1</sup>

$$y = y_0 + x - x_0 + h \cdot \tanh(k_0 \cdot (x - x_0)) \tag{1}$$

where:

$$x = \ln \frac{p}{p_0} \tag{2}$$

$$y = \ln \frac{s}{1 - s} \tag{3}$$

$$x_0 = a + b \tag{4}$$

$$h = h_0 + a \tag{5}$$

The terms  $y_0$ ,  $p_0$ ,  $k_0$ , and  $h_0$  are all constants.

The terms a and b represent the chemical and thermal affinity parameters that cause the curve to shift position.

Hence, expanded:

$$y = y_0 + x - a - b + [(h_0 + a) \cdot \tanh(k_0 \cdot (x - a - b))]$$
(6)

$$\ln \frac{s}{1-s} = y_0 + \ln \left(\frac{p}{p_0}\right) - a - b + \left[ (h_0 + a) \cdot \tanh(k_0 \cdot (\ln \left(\frac{p}{p_0}\right) - a - b)) \right] \tag{7}$$

## **Problems**

- 1. Find a when all other terms are known.
- 2. Find p when all other terms are known.

No explicit algebraic solution exist to express the equation in terms of a or p, and hence an iterative approach is required. The literature suggests using a Newton-Raphson method.

<sup>&</sup>lt;sup>1</sup>It's actually a little more complicated, as p is the combined partial pressure of oxygen and carbon monoxide, and s is the combined saturation of oxygen and carbon monoxide, but as carbon monoxide tends to zero, we can assume  $p \approx pO_2$  and  $s \approx sO_2$ .

## Find a

This method is used to find a, the chemical affinity parameter that influences the position of the curve, from a single paired measurement of p and s (e.g. from an ABG sample).

- 1. Measure p and s.
- 2. Estimate a temporary value for a. Sensible initial guesses include a=0 for a 'standard' curve ( $pH=7.40, pCO_2=5.33, cDPG=5.00$ , with no dyshaemoglobins), or a=-0.00475 for a standard curve with trace dyshaemoglobins.
- 3. Use this temporary a and the measured p to calculate a temporary s (Eqns (3) and (1)).
- 4. The difference between the temporary and the measured s allows the calculation of a new temporary parameter a, using a Newton-Raphson iteration procedure.
- 5. Iteration continues until the temporary and the measured s matches within a given limit ( $\epsilon = 0.000001$ ).

## Find p

This method is used to find  $p_{50}$  by defining s=0.5, when a and b are known – either having been measured directly or calculated previously (e.g. in the above procedure)

- 1. Measure/define s, and calculate the affinity parameters a and b.
- 2. Guess a temporary p (e.g. p = 7) and calculate a temporary s (Eqns (1) and then rearrange (3)).
- 3. The difference between the temporary s and the measured/defined s allows the calculation of a new temporary p using a fast Newton-Raphson procedure.
- 4. The procedure is iterated until the difference between the temporary s and the measured/defined s is less than a given limit ( $\epsilon=0.000001$ ).